

INFECTION PREVENTION

Make Surfaces Safe.



SoChlor[®] DST

WHAT IS DST ?

DST or Detergent Sanitiser Tablets contain NaDCC together with a low foaming anionic surfactant. A solution made from a SoChlor Detergent Sanitiser Tablet will disinfect, degrease and clean at the same time offering greater efficiency of use for all institutional domestic staff.

As well as detergent builders and low-foaming surfactants for superior cleaning results, solutions prepared from SoChlor Detergent Sanitising Tablets are fast acting and have a complete spectrum of biocidal activity. Bacteria, bacterial spores, algae, fungi, protozoa and viruses are sensitive to their effects.

Effective on Biofilms: SoChlor DST's dual action is also proven to be one of the most effective products on hard surface biofilms.

Certified with test standards: EN 1276 & EN 13727 Bactericidal, EN 13697 Bactericidal & Fungicidal, EN 1650 Fungicidal, EN 14348 Mycobacteria, EN 13704 Sporocidal & EN 14476 Virucidal.

COVID-19 - DEEP CLEAN

SoChlor[®] DST

- ✓ Dual-action / Clean & Disinfect. No pre-clean required.
- ✓ NaDCC recommended by CDC, WHO & UK Department of Health for effective disinfection of Coronavirus.
- ✓ EN 14476 certified.
- ✓ 1,000ppm av.cl recommended by Public Health England



Fully Effective Against Coronavirus

www.gvhealth.com

DST OFFERS:

- Broad spectrum efficacy suitable for a wide variety of applications, including blood spills and sporicidal disinfection
- Easy to use formats with detailed training materials and guidance available
- Single-step cleaning and disinfection with DST



HEALTHIER SAFER HUMAN ENVIRONMENTS

GV Health Ltd, Hall House, 2 Arlington Court, Whittle Way,
Arlington Business Park, Stevenage, Hertfordshire, SG1 2FS, UK
T +44 (0)1920 463 098 F +44 (0)1920 484 664 E info@gvhealth.com
E export@gvhealth.com (outside UK)

www.gvhealth.com

HOW TO MAKE SOLUTION

Tablet format allows for exact dilution. For example; to create a 1,000ppm solution add 1 tablet to 1 litre of water.



Application

Healthcare Clean Disinfect
Outbreak Disinfection
Biohazard Spill Disinfection



Chlorine concentration parts per million (ppm)

1,000
5,000
10,000

Preparation

See tablet below. Add tablet(s) as per requirement to 1 litre of water to create chlorine dual-action solution.

Tablet Dilution per 1 litre

1
5
10

*Public Health England recommend 1,000ppm for disinfection of COVID-19

Surfaces

Wipe down all bed frames, mattress covers and bed side surfaces. Apply the solution to surface and leave to dry.

APPLICATIONS

Healthcare Infection Control


- Daily clean
- Terminal clean
- Isolation clean
- Deep clean
- Outbreak decontamination
- Biohazard spill decontamination

Public Health

- Environmental Disinfection / Decontamination
- Deep clean
- Outbreak decontamination
- Biohazard spill decontamination

Floors

Mop all floor areas and surfaces under beds. Apply the solution to surface and leave to dry.



INFECTION PREVENTION SUPPORT & TRAINING

SoChlor tablets are manufactured under GMP and ISO13485 standards ensuring the highest standards of quality and compliance. SoChlor is uniquely formulated to ensure market leading solutions for environmental disinfection and decontamination.

But manufacture is only half the story. Our infection control support teams and user tools ensure healthcare providers achieve the highest efficacy and standards for infection prevention. We provide 365 support, with implementation workshops and training for staff. We also offer bespoke posters and diluter labeling as per hospital infection control requirement.

As used by the NHS: SoChlor is available on the NHS Supply Chain via the codes below.



Find out more about SoChlor TAB range on our website.

Order Code	Description	NaDCC / Chlorine per tablet	Quantity
MFB256	SoChlor™ DST (Chlorine + Detergent release)	1.7g / 1.02g	100 Tablets per pot 6x pots per case
MFB294	SoChlor™ DST (Chlorine + Detergent release)	1.7g / 1.02g	200 Tablets per pot 6x pots per case
MFB257	SoChlor™ Diluter 1 Litre		
MFB258	SoChlor™ Diluter 2 Litre		
MFB259	SoChlor™ Diluter 5 Litre		
MFB693	SoChlor™ Spray Diffuser 1 Litre		

